

ASPIRATORS

Aspirators

High density polythene

- Robust design, graduated in litres.
- Integrally moulded spout threaded to accept replaceable polypropylene tap.
- Supplied with handle and leakproof polythene caps.
- Tap assembly increases width by 35mm.

Code	Capacity litres	Height mm	Width mm	Quantity per pack
P023/01	5	332	170	13
P023/02	10	378	205	9

P023/73 Spare black cap for 5 & 10L aspirators

P023/04 Spare screw on tap



Aspirator/Jerrican

High density polythene

- Integrally moulded handle
- Supplied with screw-in bung to allow use with or without screw-in tap
- Fitted with polythene leak-proof cap

Code	Capacity litres	Height mm	Width mm	Depth mm	Type
P025/01	25	460	330	220	without tap (jerrican)
P025/03	25	460	330	220	with tap (aspirator)

P025/02 Spare tap

P025/72 Spare polythene cap

P025/73 Spare bung

