



Specifications:

pH Range

pH 0 to 14 ($\text{Na}^+ < 0.1\text{N}$)

Response Time

10 seconds (95% Response)

from pH 7.00 to pH 4.01

Glass Membrane Impedance

20 Mohms at 25°C

GX-V Glass

Introducing our new pH glass formula with unparalleled performance. This new pH glass formula, GX-V, can replace all other formulas, streamlining inventories and eliminating wasted time trying to select the best pH electrode for the job. The GX-V formulation is perfect in low ionic strength solutions as well as high pH applications and is specially formulated to have lightning fast response time with superior accuracy and repeatability. The sample pH is registered in less than 10 seconds, the result is correct to within 0.1 pH units and the pH value is reproducible over the entire pH range from 0 to 14 pH. For the best pH glass formula on the market, **just think green.**