

Colorful & Convenient mini-centrifuge

Cubee

Producer: GeneReach biotechnology Corp.

Origin : Taiwan



瑞基海洋生物科技

GeneReach Biotechnology Corp

Features

- 3 in 1 Rotor, patented
 - 8 × 1.5~2 ml tube,
 - 8 × 0.5 ml tube,
 - 16 × 0.02 ml PCR tube
- Durability
 1. Pass 5-day Burn-in test
 2. Startup switch can be used more than 200,000
 3. Strengthen bearing of lid
- User friendly: operating by single hand
- Colorful & Cute

Specification

- Dimensions

18.3 × 15.1 × 15.3 cm (W×D×H)

- Power:

100~240 V, auto-switch, 50/60 Hz

- Weight

1.8 kg

- Rotor Speed

4,000 rpm

- Noise

- < 70 dB

3 in 1 rotor, Patented



8 × 1.5~2 ml tube
8 × 0.5 ml tube
16 × 0.02 ml PCR tube



瑞基海洋生物科技
GeneReach Biotechnology Corp.

Durable design



Longer bearing of lid



Thicken bearing

瑞基海洋生物科技
GeneReach Biotechnology Corp.

Convenient! Can be used by one hand



瑞基海洋生物科技
GeneReach Biotechnology Corp.

Packing



瑞基海洋生物科技
GeneReach Biotechnology Corp.



Total 5 color

1. Orange
2. Blue
3. Pink
4. Yellow
5. Green

瑞基海洋生物科技
GeneReach Biotechnology Corp.

Contact information

Contact: Eric Hou

Mobile : +886-919-265345

Tel : +886-4-2463-9869

Fax : +886-4-2463-8255

E-mail : eric@genereachbiotech.com

GeneReach Biotechnology Corp

No.19, Keyuan 2nd Rd., Central Taiwan Science Park ,
Taichung City 407, Taiwan