



A wide range of pouches featuring a self-adhesive strip which by removing the backing strip and folding over it's possible to create an integral seal that would withstand sterilisation. This pouch offers a combination of high strength and good barrier properties with a soft texture.

- Manufactured from tinted blue, multi-layer PET/PP film
- 68gsm medical grade paper
- Chevron shaped end with easy opening
- Printed with Class 1 process indicators compliant with ISO.11140-1
- Sterilisation indicators comply with BS EN ISO 11140-1:2005
- Manufactured under BS EN ISO.9001.2008 and BS EN ISO 13485:2003
- Conforms to BS EN ISO11607-1:2006 part 1 and BS EN 868-3 & 5:2009

- STOCK SIZES

Product Code	Size (mm)	Pack Qty
WPSS060100	60 x 100	1000
WPSS090140	90 x 140	1000
WPSS090200	90 x 200	1000
WPSS090230	90 x 230	1000
WPSS135255	135 x 255	1000
WPSS135360	135 x 360	1000
WPSS190330	190 x 330	1000
WPSS300450	300 x 450	1000