

QUANTOFIX® test strips

Rapid results directly at the point of interest



QUANTOFIX® test strips are semi-quantitative in nature. They meet the most important requirements for a modern quick test.

Rapid

Most of the QUANTOFIX® test strips can be carried out in 10-120 seconds. That makes the application very quick for the end user.

Convenient

All QUANTOFIX® tests are ready to use kits. They are precalibrated and contain all necessary equipment and reagents. As "labs in a pocket" they are very handy for the end user.

Precise

The color charts are adjusted and checked using certified standard solutions that are directly traceable to. The user can be sure to receive accurate readings whenever he tests.

Applications

| | |
|--------------------------|--|
| Aquaculture | QUANTOFIX® Ammonium QUANTOFIX® Carbonate Hardness QUANTOFIX® Nitrate/Nitrite QUANTOFIX® Multistick for aquarium owners |
| Electroplating | QUANTOFIX® Chromate QUANTOFIX® Copper QUANTOFIX® Cyanide QUANTOFIX® Nickel QUANTOFIX® Nitrate QUANTOFIX® Zinc |
| Metall workshops | QUANTOFIX® Nitrate QUANTOFIX® Nitrite QUANTOFIX® LubriCheck |
| Food processing industry | QUANTOFIX® Ascorbic Acid QUANTOFIX® Calcium QUANTOFIX® Carbonate Hardness QUANTOFIX® Chloride QUANTOFIX® Total acid QUANTOFIX® Total sugar QUANTOFIX® Nitrate/Nitrite QUANTOFIX® Peroxide QUANTOFIX® Phosphate QUANTOFIX® QUAT QUANTOFIX® Sulphite |

QUANTOFIX® test sticks at a glance

| Test | Gradation | REF |
|-----------------------------------|--|-------|
| QUANTOFIX® Active oxygen | 0 · 4 · 8 · 15 · 25 mg/L KMPS | 91349 |
| QUANTOFIX® Aluminium | 0 · 5 · 20 · 50 · 200 · 500 mg/L Al ₃₊ | 91307 |
| QUANTOFIX® Ammonium | 0 · 10 · 25 · 50 · 100 · 200 · 400 mg/L NH ₄₊ | 91315 |
| QUANTOFIX® Arsenic Sensitive | 0 · 0.005 · 0.01 · 0.025 · 0.05 · 0.1 · 0.25 · 0.5 mg/L As _{3+/5+} | 91345 |
| QUANTOFIX® Arsenic 10 | 0 · 0.01 · 0.025 · 0.05 · 0.1 · 0.5 mg/L As _{3+/5+} | 91334 |
| QUANTOFIX® Arsenic 50 | 0 · 0.05 · 0.1 · 0.5 · 1.0 · 1.7 · 3.0 mg/L As _{3+/5+} | 91332 |
| QUANTOFIX® Ascorbic acid | 0 · 50 · 100 · 200 · 300 · 500 · 1000 · 2000 mg/L ascorbic acid | 91314 |
| QUANTOFIX® Calcium | 0 · 10 · 25 · 50 · 100 mg/L Ca ₂₊ | 91324 |
| QUANTOFIX® Carbonate hardness | 0 · 3 · 6 · 10 · 15 · 20 °d | 91323 |
| QUANTOFIX® Chloride | 0 · 500 · 1000 · 1500 · 2000 · ≥3000 mg/L Cl ⁻ | 91321 |
| QUANTOFIX® Chlorine | 0 · 1 · 3 · 10 · 30 · 100 mg/L Cl ₂ | 91317 |
| QUANTOFIX® Chlorine Sensitive *CE | 0 · 0.1 · 0.5 · 1 · 3 · 10 mg/L Cl ₂ | 91339 |
| QUANTOFIX® Chromate | 0 · 3 · 10 · 30 · 100 mg/L CrO ₄₂₋ | 91301 |
| QUANTOFIX® Cobalt | 0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 mg/L Co ₂₊ | 91303 |
| QUANTOFIX® Copper | 0 · 10 · 30 · 100 · 300 mg/L Cu _{1+/2+} | 91304 |
| QUANTOFIX® Cyanide | 0 · 1 · 3 · 10 · 30 mg/L CN ⁻ | 91318 |
| QUANTOFIX® EDTA | 0 · 100 · 200 · 300 · 400 mg/L EDTA | 91335 |
| QUANTOFIX® Formaldehyde | 0 · 10 · 20 · 40 · 60 · 100 · 200 mg/L HCHO | 91328 |
| QUANTOFIX® Glucose | 0 · 50 · 100 · 250 · 500 · 1000 · 2000 mg/L Glucose | 91348 |
| QUANTOFIX® Glutaraldehyde CE* | 0 · 0.5 · 1 · 1.5 · 2 · 2.5 % Glutaraldehyde | 91343 |
| QUANTOFIX® LubriCheck | 0 · 15 · 50 · 75 · 130 · 200 mmol/L KOH | 91336 |
| QUANTOFIX® Molybdenum | 0 · 5 · 20 · 50 · 100 · 250 mg/L Mo ₆₊ | 91325 |
| QUANTOFIX® Nickel | 0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 mg/L Ni ₂₊ | 91305 |
| QUANTOFIX® Nitrate 100 | 0 · 5 · 10 · 25 · 50 · 75 · 100 mg/L NO ₃₋ 0 · 0.5 · 2 · 5 · 10 · 25 · 50 mg/L NO ₂₋ | 91351 |
| QUANTOFIX® Nitrate/Nitrite | 0 · 10 · 25 · 50 · 100 · 250 · 500 mg/L NO ₃₋ / 0 · 1 · 5 · 10 · 20 · 40 · 80 mg/L NO ₂₋ | 91313 |
| QUANTOFIX® Nitrite | 0 · 1 · 5 · 10 · 20 · 40 · 80 mg/L NO ₂₋ | 91311 |
| QUANTOFIX® Nitrite 3000 | 0 · 100 · 300 · 600 · 1000 · 2000 · 3000 g/L NO ₂₋ | 91322 |
| QUANTOFIX® Nitrite/pH | 0 · 1 · 5 · 10 · 20 · 40 · 80 mg/L NO ₂ pH 6.0 · 6.4 · 6.7 · 7.0 · 7.3 · 7.6 · 7.9 · 8.2 · 8.4 · 8.6 · 8.8 · 9.0 · 9.3 · 9.6 | 91338 |

| | | |
|--|---|-------|
| QUANTOFIX® Peracetic acid 50 CE* | 0 · 5 · 10 · 20 · 30 · 50 mg/L PAA | 91340 |
| QUANTOFIX® Peracetic acid 500 CE* | 0 · 50 · 100 · 200 · 300 · 400 · 500 mg/L PAA | 91341 |
| QUANTOFIX® Peracetic acid 2000 CE* | 0 · 500 · 1000 · 1500 · 2000 mg/L PAA | 91342 |
| QUANTOFIX® Peroxide 25 | 0 · 0,5 · 2 · 5 · 10 · 25 mg/L H ₂ O ₂ | 91319 |
| QUANTOFIX® Peroxide 100 CE* | 0 · 1 · 3 · 10 · 30 · 100 mg/L H ₂ O ₂ | 91312 |
| QUANTOFIX® Peroxide 1000 | 0 · 50 · 150 · 300 · 500 · 800 · 1000 mg/L H ₂ O ₂ | 91333 |
| QUANTOFIX® Phosphate | 0 · 3 · 10 · 25 · 50 · 100 mg/L PO ₄₃₋ | 91320 |
| QUANTOFIX® Potassium | 0 · 200 · 400 · 700 · 1000 · 1500 mg/L K ⁺ | 91316 |
| QUANTOFIX® QUAT | 0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 mg/L benzalkoniumchloride | 91337 |
| QUANTOFIX® Silver | 1 · 2 · 3 · 5 · 7 · 10 g/L Ag ⁺ | 91350 |
| QUANTOFIX® Total acid | 0 · 2.0 · 2.5 · 3.0 · 3.5 · 4.0 · 5.0 g/L citric acid | 91353 |
| QUANTOFIX® Total Iron 1000 | 0 · 5 · 20 · 50 · 100 · 250 · 500 · 1000 mg/L Fe _{2+/3+} | 91330 |
| QUANTOFIX® Total Iron 100 | 0 · 2 · 5 · 10 · 25 · 50 · 100 mg/L Fe _{2+/3+} | 91344 |
| QUANTOFIX® Total sugar | 0 · 55 · 100 · 250 · 400 · 600 · 800 mg/L Fructose / Glucose | 91352 |
| QUANTOFIX® Sulfate | <200 · >400 · >800 · >1200 · >1600 mg/L SO ₄₂₋ | 91329 |
| QUANTOFIX® Sulfite | 0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 mg/L SO ₃₂₋ | 91306 |
| QUANTOFIX® Tin | 0 · 10 · 25 · 50 · 100 · 250 · 500 mg/L Sn ₂₊ | 91309 |
| QUANTOFIX® Zinc | 0 · 2 · 5 · 10 · 25 · 50 · 100 mg/L Zn ₂₊ | 91310 |
| QUANTOFIX® Multisticks for aquarium owners | 0 · 5 · 10 · 15 · 20 · 25 °d total hardness 0 · 3 · 6 · 10 · 15 · 20 °d carbonate hardness pH 6.4 · 6.8 · 7.2 · 7.6 · 8.0 · 8.4 | |