

## Sterilisation Reels



- Suitable for Steam Sterilisation 
- Suitable for ETO Sterilisation 
- Suitable for Formaldehyde 
- CE Marked 
- BSI Kite Mark 
- Single Use 

These have been designed for packaging extra long products or creating individual pouches of special lengths.

- Medical grade paper 60g
- Tinted, multi-layer, heat-sealable PET/PP film
- Available in flat or gusseted for larger objects
- Transparent and tinted blue as standard
- Printed with Class 1 process indicators compliant with ISO.11140-1
- Sterilisation indicators comply with BS EN ISO 11140-1:2005
- Manufactured under BS EN ISO.9001.2008 and BS EN ISO 13485:2003
- Conforms to BS EN ISO11607-1:2006 part 1 and BS EN 868-3 & 5:2009
- Shelf life 5 years from manufacture date
- Storage in cool and dark place

### Flat Stock Sizes (200m length)

Product Code	Size (mm)	Pack Qty
WR050	50	1
WR075	75	1
WR100	100	1
WR150	150	1
WR200	200	1
WR250	250	1
WR300	300	1
WR400	400	1