

Specialty Products for Color Coding

Shamrock's Multi-Purpose Labeling Tape

Provides the Perfect Solution to a Wide Variety of Labeling Needs.

- Versatile
- Dependable
- Transferable
- Acid Resistant
- Waterproof
- Oilproof



IDEAL FOR COLOR CODING

Adheres to any dry surface – glass, plastic, metal, paper and rubber.

Sticks throughout the most demanding procedures – is oilproof, waterproof and acid resistant – can be steam autoclaved.

Holds in a wide range of temperatures – transfers cleanly from one surface to another, leaving no messy residue.

Accepts handwritten notations with pen, pencil or marker, or order pre-printed to your specifications.

31 Colors

1 White	11 Rose	22 Maroon
2 Yellow	12 Violet	23 Dark Brown
3 Green	13 Lavender	24 Dark Blue
4 Red	14 Gray	25 Dark Green
5 Orange	15 Lime	26 Black
6 Blue	16 Aqua	27 Fl. Red
7 Pink	17 Chartreuse	28 Fl. Orange
8 Tan	18 Gold	29 Fl. Green
9 Salmon	19 Silver	30 Fl. Chartreuse
10 Copper	20 Magenta	31 Fl. Pink
	21 Dark Red	



Individual Rolls of Shamrock Labeling Tape are available in your choice of:

- 31 Colors
- 2 Lengths
- 5 Widths

Width	500" Cat. No.	2160" Cat. No.
1/2"	ST-12	ST-120
3/4"	ST-34	ST-340
1"	ST-10	ST-100
1 1/2"	ST-11	ST-110
2"	ST-20	ST-200

Rainbow Packs of tape rolls in white, yellow, green, red, orange, blue and pink are also available.

- 7 Colors
- 500" Tape Rolls
- Your Choice of 3 Widths

Cat. No.	Width	Rolls per Pack
ST-12-RNBW	1/2"	24
ST-34-RNBW	3/4"	16
ST-10-RNBW	1"	12



MULTIPLE ROLL TAPE DISPENSERS

Easy to fill dispensers, available in four sizes, offer convenient access to tape from a central location.



HAND DISPENSERS

Durable plastic, can be reused many times.

