



## Hydrion S/R Dispenser 1.0-14.0

Catalog#:94

Hydrion Insta-Chek pH Paper, the ultimate in wide range, general purpose pH papers has a distinct color match at each full pH

unit:[1.0][2.0][3.0][4.0][5.0][6.0][7.0][8.0][9.0][10.0][11.0][12.0][13.0][14.0]

Hydrion Spectral Single Roll pH Paper provides:

A direct reading pH paper with the widest range available in a single paper

A simple, single color match for all pH values

Distinct and easily matched color separations

Accurate, reproducible readings with an immediate response

Hydrion Spectral Single Roll pH Paper is a general use pH paper. A sample of its diverse pH applications include:

Classroom Demonstrations

Hospital Laboratories

General Laboratories

Industrial Testing

Food Service

Safety Testing

Emergency Testing

Agricultural Testing

Each carton contains 10 kits, each consisting of a 15-foot roll of test paper and matching color chart. Provides for approximately 1000 tests.