

Item	Packaging(pcs/pk)		Description	
M047PTFE022H		Pore:0.22μm	Diameter:47mm	200
M047PTFE045H		Pore:0.45μm	Diameter:47mm	200
M047PTFE010H		Pore:1.00μm	Diameter:47mm	200
M047PTFE30H		Pore:3.0μm	Diameter:47mm	200
M047PTFE500H		Pore:5.0μm	Diameter:47mm	200

Hydrophobic PTFE Membrane Filter

PTFE membrane with supporting layer polyester or polypropylene

Suitable for applications involving aggressive organic solvents, strong acids, and alkalis

Pore size: 0.22μm-5μm, Diameter: 47mm

Hydrophobic nature of the membrane has applications for air and gas sterilization

High temperature resistance

