




DURAN® GL 45 Bottle 2000ml with GL 45 Side Neck





DURAN® GL 45 Bottle 2000ml with GL 45 Side Neck



APPLICATION

- ▶ HPLC Reservoir Bottle
- ▶ Instrument Reagent Supply Bottle

FEATURES & BENEFITS

- Convenient – Allow filling without main cap removal
- Manufactured from the material of choice for HPLC solvent bottles: borosilicate 3.3 glass offering a very low inorganic ion extractables profile
- GL 45 thread fits standard HPLC connection caps
- Useful 2 Liter Capacity
- Retrace Code, and certificate for manufacturing lot traceability
- Durable white enamel volume graduations



What is HPLC “Purging”?

When setting up an HPLC system, the aim of the purge is simply to flush through all the lines so that any remaining solvent in them from a previous analysis or wash is replaced with the new mobile phase.

