


01 LABORATORY GLASS BOTTLES AND ACCESSORIES

DURAN® Screw Cap with aperture  
from PBT, red

Suitable for silicone seal for piercing (Septa). More chemically resistant than PP cap.

Typical applications: Injection or removal of media.



A  
121 °C

Tmax.  
180 °C

| Cat. No.  | DIN Thread (GL) | Aperture bore d (OD) (mm) | d (OD) (mm) | h (mm) | Pack Unit |
|-----------|-----------------|---------------------------|-------------|--------|-----------|
| 292270508 | 14              | 9.5                       | 20          | 17     | 10        |
| 292270602 | 18              | 11                        | 23          | 20     | 10        |
| 292270902 | 25              | 15                        | 33          | 23     | 10        |
| 292270808 | 32              | 20                        | 42          | 26     | 10        |
| 292271007 | 45              | 34                        | 54          | 28     | 10        |