



Pyrex® Centrifuge tube pear shape, ASTM

ASTM Pear shape centrifuge tubes

Complying to ASTM D96 and D-1966, these 100ml pear shape centrifuge tubes are suitable for the determination of water and sediment in petroleum products.

- Complies with ASTM D96 and D-1966
- For the determination of water and sediment in petroleum products
- Graduated in white enamel: - Lower stem graduated 0 to 1.5ml in 0.1ml divisions - Upper body graduated 1.5 to 5ml in 0.5ml divisions and from 5 to 10ml in 1ml



Product Specification

Item Code	Approx overall length mm	Capacity ml	Neck ID mm	Pack quantity
3440/02M	157	100	17	2