

## DURAN® MULTIPOINT CONNECTOR CAPS GL 25

2, 3 or 4 port connector caps for GL 25 DURAN® bottles

### Applications

- Low volume sampling
- Media supply for small bioreactors
- Fed-batch bioreactor
- Perfusion circuit for cell culture
- Oxygenation / aeration of very small samples
- Use in a chemistry micro-reactor

**NEW**



### Features

- Caps have a GL 25 thread that fit the smaller 10 or 25 ml DURAN® bottles
- Barbed male hose connector are suitable for tubing with a 3.2 mm (1/8") inner diameter (such as PharMed® BPT Tubing\*)
- Available with 2, 3, or 4 surgical grade stainless steel (316L) barbed tubing ports
- Construction materials are certified as approved for food contact
- High temperature resistance (up to 180 °C)
- Autoclavable (@ 121 °C / 15 minutes)
- Reusable with a long service life

\*PharMed® is a registered trademark of Saint-Gobain

## Ordering information

Cat. No.	Description	Pack Quantity
29 260 13 01 <b>NEW</b>	DURAN® 2-port Connector caps GL 25, red PBT screw cap, PTFE insert and 2 stainless steel barbs	1
29 261 13 02 <b>NEW</b>	DURAN® 3-port Connector caps GL 25, red PBT screw cap, PTFE insert and 3 stainless steel barbs	1
29 262 13 03 <b>NEW</b>	DURAN® 4-port Connector caps GL 25, red PBT screw cap, PTFE insert and 4 stainless steel barbs	1

## Recommended compatible DURAN® bottles

Cat. No.	Description	Capacity ml	Pack Quantity
21 801 08 02	DURAN® Clear GL 25 bottle supplied without screw cap. Retrace coded.	10 ml	10
21 806 08 07	DURAN® Amber GL 25 bottle supplied without screw cap. Retrace coded.	10 ml	10
21 805 08 06	DURAN® Clear protect safety coated GL 25 bottle supplied without screw cap. Retrace coded.	10 ml	10
21 801 14 04	DURAN® Clear GL 25 bottle supplied without screw cap. Retrace coded.	25 ml	10
21 806 14 09	DURAN® Amber GL 25 bottle supplied without screw cap. Retrace coded.	25 ml	10
10 926 76	DURAN® Clear protect safety coated GL 25 bottle supplied without screw cap. Retrace coded.	25 ml	10
21 806 14 33	DURAN® Amber protect safety coated GL 25 bottle supplied without screw cap. Retrace coded.	25 ml	10

These products are available from your specialist laboratory dealer.