

The 80 mm wide outer diameter of the bottle neck permits easy filling and pouring out of powders and viscous substances. With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Available as sole bottle and completed bottle with blue quick release closure (PP¹, integral lip seal) and pouring ring (PP¹) for drip-free pouring and clean, safe working. Service temperature limit of closure and pouring ring: +140 °C. Special thread means opening takes less than a turn.


Typical applications: storage, transport, safekeeping and sampling of substances, easy to use with granulated material, powders and viscous media, sampling of hot media.

Cat. No. with cap	Cat No. without cap	Capacity ml	Thread GLS	d mm	h mm	Pack/ Quantity
21 860 36 5	21 860 36	250	80	95	105	10
11 126 27	11 783 92	500	80	101	152	10
11 127 13	11 784 24	1 000	80	101	222	10
11 127 15	11 784 25	2 000	80	136	252	10
21 860 69 5	21 860 69	3 500	80	160	275	1
11 139 49	11 784 26	5 000	80	182	314	1
11 139 50	11 784 27	10 000	80	227	389	1
11 139 51	11 784 28	20 000	80	288	484	1
12 002 65 ²		30 000	80	340	548	1
12 001 54 ²		50 000	80	400	590	1

¹ For chemical and thermal resistance, see page 49.

² Laboratory bottles 30 and 50 l available upon request. Please contact us.

DURAN® GLS 80® LABORATORY BOTTLE, WIDE NECK with GLS 80® thread



The 80 mm wide outer diameter of the bottle neck permits easy filling and pouring out of powders and viscous substances. With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Available as sole bottle and completed bottle with blue quick release closure (PP¹, integral lip seal) and pouring ring (PP¹) for drip-free pouring and clean, safe working. Service temperature limit of closure and pouring ring: +140 °C. Alongside easy handling, UV protection up to 500 nm. Unchanged DURAN® properties within the bottle, as colouration is only to the outer surface. Very uniform, durable and chemically resistant amber colour due to use of innovative technology.

Typical applications: storage, transport and safekeeping of light-sensitive substances, easy to use with granulated material, powders and viscous media.

Cat. No. with cap	Cat No. without cap	Capacity ml	Thread GLS	d mm	h mm	Pack/ Quantity
21 866 36 5	21 866 36	250	80	95	105	10
11 601 46	11 784 29	500	80	101	152	10
11 601 47	11 784 30	1 000	80	101	222	10
11 601 48	11 784 31	2 000	80	136	252	10
21 866 69 5	21 866 69	3 500	80	160	275	1
11 601 49	11 784 32	5 000	80	182	314	1
11 601 50	11 784 33	10 000	80	227	389	1
11 601 51	11 784 34	20 000	80	288	484	1

¹ For chemical and thermal resistance, see page 49.

DURAN® GLS 80® LABORATORY BOTTLE, WIDE NECK, AMBER with GLS 80® thread

