



YOUTILITY – TOUCH THE FUTURE

THE NEW DURAN® LABORATORY BOTTLE SYSTEM

Complete DURAN® YOUTILITY bottle sets GL 45 and DURAN® YOUTILITY Screw Caps, Pouring Rings, Labels and Bottle Tag selection:

DURAN® YOUTILITY bottle GL 45, clear glass, complete				
Cat. No.	Nominal Volume (ml)	Height (mm)	Outer diameter (mm)	Pack size
21 881 28 54	125	124	55	4
21 881 36 53	250	158	66	4
21 881 44 52	500	193	78	4
21 881 54 57	1000	253	93	4
DURAN® YOUTILITY bottle GL 45, amber glass, complete				
21 886 28 59	125	124	55	4
21 886 36 58	250	158	66	4
21 886 44 57	500	193	78	4
21 886 54 53	1000	253	93	4

DURAN® YOUTILITY accessories and replacement parts			
Cat. No.	Description	Material and Colour	Pack size
29 229 28 02	DURAN® YOUTILITY Screw Cap GL 45	Food grade Polypropylene (PP), cyan	10
29 241 28 08	DURAN® YOUTILITY Pouring Ring GL 45	Polypropylene (PP), cyan	16
29 243 29 04	DURAN® Bottle Tags GL 45	Silicone, eight colours. Two Bottle Tags are supplied in each colour	16
29 401 00 06	DURAN® YOUTILITY 200 self-adhesive label sets	White polyester labels in a dispenser box	1 x 200 sets

For more information about the DURAN® YOUTILITY range get in touch at: duran-youtility.com