


02 BOILING FLASKS AND GENERAL LABORATORY GLASSWARE

DURAN® Beaker

high form, with spout

With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. With spout for clean pouring. Uniform wall thickness distribution makes these beakers ideal for heating applications.



Cat. No.	Capacity (mL)	d (OD) (mm)	h (mm)	Pack Unit
211161704	50	38	70	10
211162409	100	48	80	10
211162906	150	54	95	10
211163602	250	60	120	10
211164101	400	70	130	10
211164804	600	80	150	10
211165303	800	90	175	10
211165406	1 000	95	180	10
211166308	2 000	120	240	10
211166805	3 000	135	280	2