

With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. Spout for clean pouring. Uniform wall thickness distribution makes these beakers ideal for heating applications.

DURAN® Beaker

low form, with spout

Cat. No.	Capacity (mL)	d (OD) (mm)	h (mm)	Remark	Pack Unit
211060701	5	22	30	Without graduation.Without Retrace Code.	10
211060804	10	26	35	Without graduation.Without Retrace Code.	10
211061406	25	34	50		10
211061706	50	42	60		10
211062402	100	50	70		10
211062908	150	60	80		10
211063604	250	70	95		10
211064103	400	80	110		10
211064806	600	90	125		10
211065305	800	100	135		10
211065408	1 000	105	145		10
211066301	2 000	132	185		10
211066807	3 000	152	210		4
211067306	5 000	170	270		1
211068602	10 000	217	350	Non-DIN/ISO size.	1

