

Funnels, filter, vacuum Büchner



| Büchner vacuum filtration funnel with glass sintered disc | Manufactured from borosilicate glass for the high resistance to chemical attack | Useful where conventional papers would be destroyed by acids and alkalis | With a ° ground glass cone to allow direct filtration into Büchner flasks

Nom. capacity ml	Porosity grade	Disc diam. mm	Diam. rim mm	of Cone size
30	3	30	40	19/26
30	4	30	40	19/26
80	1	40	50	24/29
80	2	40	50	24/29
80	3	40	50	24/29
80	4	40	50	24/29
125	1	60	71	24/29
125	2	60	71	24/29
125	3	60	71	24/29
125	4	60	71	24/29
500	2	95	106	29/32
500	3	95	106	29/32
500	4	95	106	29/32