

DURAN® Crystallizing Dish

with and without spout



DIN
12337

DIN
12338

A
121 °C

Cat. No.	Capacity (mL)	d (OD) (mm)	h (mm)	Pack Unit
with spout, DIN 12 338				
213112401	20	40	25	10
213113209	40	50	30	10
213113406	60	60	35	10
213113809	100	70	40	10
213114102	150	80	45	10
213114402	300	95	55	10
213114908	500	115	65	10
213115407	900	140	75	10
213115904	2 000	190	90	10
213116309	3 500	230	100	10
without spout, DIN 12 337				
213132403	20	40	25	10
213133202	40	50	30	10
213133408	60	60	35	10
213133802	100	70	40	10
213134104	150	80	45	10
213134404	300	95	55	10
213134901	500	115	65	10
213135409	900	140	75	10
213135906	2 000	190	90	10
213136302	3 500	230	100	10