




DURAN® MEASURING CYLINDER, CLASS A

- Blue graduation, with main points ring graduation and hexagonal base
- Certificate of conformity "H"
- [With batch certificate](#)



Ordering information:

Cat. No.	Capacity ml	Accuracy limit ± ml	Graduation divisions ml	h mm	d mm	Pack Quantity
21 390 07 01	5	0,05	0,1	112	12	2
21 390 08 04	10	0,1	0,2	137	14	2
21 390 14 06	25	0,25	0,5	167	21	2
21 390 17 06	50	0,5	1	196	25	2
21 390 24 02	100	0,5	1	256	29	2
21 390 36 04	250	1	2	331	39	2
21 390 44 03	500	2,5	5	360	53	2
21 390 54 08	1000	5	10	460	65	1
21 390 63 01	2000	10	20	500	85	1