


Best.-Nr. Cat. No.	DN	h mm	ID Flansch ID Flange mm	AD Flansch OD Flange mm	Verpackungseinheit/ Stück Pack/Quantity
24 410 46	100	75	119	153 +/- 2	1
24 410 57	150	98	172	215 +/- 2	1
24 410 61	200	107	224	270 +/- 2	1
24 410 66	250	122	274	320 +/- 2	1
24 410 69	300	150	332	380 +/- 2	1

**DURAN® EXSIKKATOR-DECKEL
mit Knopf, passend für alle Unterteil-Typen**

*DURAN® DESICCATOR LID
with knob, suitable for all types of bases*




DIN
12491

Best.-Nr. Cat. No.	DN	h mm	ID Flansch ID Flange mm	AD Flansch OD Flange mm	NS	Verpackungseinheit/ Stück Pack/Quantity
24 420 46	100	62	119	153 +/- 2	24/29	1
24 420 57	150	85	172	215 +/- 2	24/29	1
24 420 61	200	94	224	270 +/- 2	24/29	1
24 420 66	250	109	274	320 +/- 2	24/29	1
24 420 69	300	137	332	380 +/- 2	24/29	1

**DURAN® EXSIKKATOR-DECKEL
für normale NS-Hähne (NS 24/29),
Typ NOVUS, passend für alle Unterteil-Typen**

*DURAN® DESICCATOR LID
for standard ground joint stopcocks (NS 24/29),
type NOVUS, suitable for all types of bases*



DIN
12491