


Best.-Nr. Cat. No.	DN	h mm	h <sub>1</sub> mm	ID Flansch ID Flange mm	AD Flansch OD Flange mm	ca. Volumen approx. Vol. l	Verpackungseinheit/ Stück Pack/Quantity
24 770 46	100	112	58	119	153 +/- 2	0,7	1
24 770 57	150	154	81	172	215 +/- 2	2,4	1
24 770 61	200	202	115	224	270 +/- 2	5,8	1
24 770 66	250	235	120	274	320 +/- 2	10,5	1
24 770 69	300	283	150	332	380 +/- 2	18,5	1

**DURAN® EXSIKKATOR-UNTERTEIL**  
mit Planflansch, ohne Anschluss,  
passend für alle Deckel-Typen

*DURAN® DESICCATOR BASE*  
with flat flange, no outlet,  
suitable for all types of lids



DIN  
12491

Best.-Nr. Cat. No.	DN	h mm	h <sub>1</sub> mm	ID Flansch ID Flange mm	AD Flansch OD Flange mm	ca. Volumen approx. Vol. l	Verpackungseinheit/ Stück Pack/Quantity
24 773 61	200	202	115	224	270 +/- 2	5,8	1

**DURAN® EXSIKKATOR-UNTERTEIL**  
mit Ringnutflansch,  
passend für alle Deckel-Typen

*DURAN® DESICCATOR BASE*  
with ring-grooved flange,  
suitable for all types of lids

