


DURAN® Desiccator Lid

with knob, suitable for all types of bases



Cat. No.	DN	ID Flange (mm)	Outer diameter (AD) Flange (mm)	h (mm)	Pack Unit
244104607	100	119	153 ± 2	75	1
244105706	150	172	215 ± 2	98	1
244106102	200	224	270 ± 2	107	1
244106608	250	274	320 ± 2	122	1
244106908	300	332	380 ± 2	150	1

סי.י.מדה בע"מ

יבוא ושיווק ציוד ומכשור למעבדות
WWW.SIMADA.CO.IL 03-5591008

