

Cat. No.	DN	ID Flange (mm)	Outer diameter (AD) Flange (mm)	h (mm)	Neck	Pack Unit
244204605	100	119	153 ± 2	62	24/29	1
244205704	150	172	215 ± 2	85	24/29	1
244206109	200	224	270 ± 2	94	24/29	1
244206606	250	274	320 ± 2	109	24/29	1
244206906	300	332	380 ± 2	137	24/29	1

ס.י.מדה בע"מ


יבוא ושיווק ציוד ומכשור למעבדות

WWW.SIMADA.CO.IL 03-5591008



DURAN® Desiccator Lid

for standard ground joint stopcocks
(NS 24/29), type NOVUS, suitable for all
types of bases



DIN ISO
13130