



Adapters, splash head, rotary evaporator, anti-climb



- For use with rotary evaporators
- Protects against liquid ingress into the condenser from 'bumping' solvents
- Useful for liquids that have tendency to foam
- Splash heads feature a 29/32 socket

Product Code	Cone size	Pack qty.
SH15/19	19/26	1
SH15/24	24/29	1
SH15/29	29/32	1



Adapters, splash head, rotary evaporator, with sinter



- For use with rotary evaporators
- Prevents possible contamination of condenser
- Features a coarse porosity 1 sinter to reduce carry over
- Splash heads feature a 29/32 socket

Product Code	Cone size	Pack qty.
SH16/19	19/26	1
SH16/24	24/29	1
SH16/29	29/32	1



Adapters, steam distillation head



- A splash-head with integral stem for the introduction of steam into flasks
- Stem is connected to steam supply via detachable screwthread safety connector

Product Code	Cone size to fit flask	Cone size to fit condenser	Pack qty.
SH7/23/SC	24/29	19/26	1



Screwthread connector, cat no: 4510/02, accepts 9mm bore flexible tubing. Spare screwthread connectors are available on page 96.

