

Splash heads. PRECISO.

Vertical cone ST	Side cone ST	Vertical socket ST	
------------------	--------------	--------------------	--

Straight

14/23	—	14/23	1.420.000
19/26	—	19/26	1.420.001
24/29	—	19/26	1.420.002
24/29	—	24/29	1.420.003
29/32	—	29/32	1.420.004


Bent

14/23	14/23	—	1.425.001
19/26	14/23	—	1.425.002
19/26	19/26	—	1.425.003
24/29	19/26	—	1.425.004
24/29	24/29	—	1.425.005
29/32	19/26	—	1.425.006
29/32	29/32	—	1.425.007


Distilling adapter.

150 ml, DURAN, PRECISO. With cone and socket ST 29/32, to be placed onto rotary evaporators.

1.430.150

 n
neuw

 n
neuw

Rotary evaporator safety trap, anti-foaming trap. DURAN, PRECISO. With cone and socket ST 29/32, to be placed between evaporating flask and condenser for use with rotary evaporators. Use where foaming is likely to occur. Relieves superheating and bumping problems. Interchangeable ST 29/32 ground joint.

1.430.200


Suction tubes. PRECISO.

Socket ST	Straight	Bent	With ST-stop cock straight
14/23	1.662.014	1.664.014	1.666.014
19/26	1.662.019	1.664.019	1.666.019
24/29	1.662.024	1.664.024	1.666.024
29/32	1.662.029	1.664.029	1.666.029


Suction tubes. PRECISO.

ST	Straight	Bent
Without stopcock		
14/23	1.680.014	1.681.014
19/26	1.680.019	1.681.019
24/29	1.680.024	1.681.024
29/32	1.680.029	1.681.029

With stopcock

14/23	1.683.014	1.682.014
19/26	1.683.019	1.682.019
24/29	1.683.024	1.682.024
29/32	1.683.029	1.682.029

With needle valve stopcock

14/23	—	1.684.014
19/26	—	1.684.019
24/29	—	1.684.024
29/32	—	1.684.029

 n
neuw

Instead of the conventional hose connection made of glass we offer our SVS-hose connection cat. no. 0.185.001. Surcharge see price list. (Please add *.***.***S behind the article no.).

Instead of the conventional hose connection made of glass we offer the components with thread RD 14. Surcharge see price list. (Please add *.***.***Rd behind the article no.).

