

05 GLASS FILTRATION APPARATUS AND ACCESSORIES

Heavy walled for vacuum use. Does not conform to the "equipment and product safety regulations".

Typical application: separations by vacuum filtration.

Cat. No.	Capacity (mL)	d (OD) (mm)	d ₁ (OD) (mm)	h (mm)	Hose connection d (OD) (mm)	Pack Unit
211916802	3 000	170	58	295	11	1
211917301	5 000	185	68	360	11	1
211918606	10 000	240	70	420	11	1
211918803	15 000	255	70	500	11	1
211919105	20 000	290	70	535	11	1

DURAN® Filtering Flask with glass hose connection

bottle shape

