


## DURAN® Filtering Flask with glass hose connection

Erlenmeyer shape



Due to the heavy wall thickness the apparatus is vacuum-tight. Does not conform to the "equipment and product safety regulations".

Typical application: separations by vacuum filtration.

Cat. No.	Capacity (mL)	d (OD) (mm)	d <sub>1</sub> (OD) (mm)	h (mm)	Hose connection d (OD) (mm)	Pack Unit
212012409	100	64	24	105	11	10
212013602	250	85	34	155	11	10
212014401	500	105	34	185	11	10
212015406	1 000	135	45	230	11	10
212016308	2 000	166	60	255	11	1