

02 BOILING FLASKS AND GENERAL LABORATORY GLASSWARE

Uniform wall thickness distribution makes these flasks ideal for heating applications. Flat base means flasks can be set down without a supporting ring. Flasks with a neck diameter of 65 mm have a reinforced rim.

DURAN® Flat Bottom Flask Narrow Neck

with beaded rim

Cat. No.	Capacity (mL)	d (OD) (mm)	d ₁ (OD) (mm)	h (mm)	Remark	Pack Unit
217111708	50	51	26	90		10
217112404	100	64	26	105		10
217113606	250	85	34	138		10
217114405	500	105	34	163		10
217115401	1 000	131	42	190		10
217116406	2 000	166	42	250	Non-DIN ISO size.	10
217116809	3 000	185	50	250	Non-DIN ISO size.	1
217117102	4 000	207	50	275		1
217117308	5 000	223	50	290	Non-DIN ISO size.	1
217117608	6 000	237	65	315		1
217118604	10 000	280	65	360		1

