


02 BOILING FLASKS AND GENERAL LABORATORY GLASSWARE

DURAN® Round Bottom Flask Narrow Neck

with beaded rim



DIN ISO
1773

A
121 °C

USP
Standard

Uniform wall thickness distribution makes these flasks ideal for heating applications. The geometry permits very uniform heating. Flasks with a neck diameter of 65 mm or more have a reinforced rim.

Cat. No.	Capacity (mL)	d (OD) (mm)	d ₁ (OD) (mm)	h (mm)	Remark	Pack Unit
217211706	50	51	26	95		10
217212402	100	64	26	110		10
217213604	250	85	34	144		10
217214403	500	105	34	168		10
217215408	1 000	131	42	200		10
217216404	2 000	166	42	250	Non-DIN ISO size.	10
217216807	3 000	185	50	260	Non-DIN ISO size.	1
217217109	4 000	207	52	290		1
217217306	5 000	223	50	305	Non-DIN ISO size.	1
217217709	6 000	236	51	355	Non-DIN ISO size.	1
217218602	10 000	279	65	380		1
217218705	12 000	295	65	380	Non-DIN ISO size. Conforms to ASTM E 1403.	1
217219101	20 000	345	76	515	Conforms to ASTM E 1403.	1