



filter funnel

Good thermal shock and chemical resistance. Filter funnels mate to the filtering flask via a conical rubber seal (GUKO).

Typical applications:

Qualitative inorganic analysis and preparative chemistry.

Cat. No.	Porosity	Capacity ml	Disk d mm	d ₁ mm	h mm	ID mm
25 852 0X *	1-5	50	35	10	130	40
25 852 1X *	1-5	75	45	10	130	56
25 852 2X *	1-5	125	60	10	140	71
25 852 3X *	1-5	500	95	22	200	106
25 852 4X *	1-5	1 000	120	22	320	135
25 852 6X *	1-4	4 000	175	30	400	202

יש לציין צפיפות הסינטר לפי הטבלה:

Porosity	New identification mark ISO 4793	Nominal max.pore size (μm)
0	P 250	160-250
1	P 160	100-160
2	P 100	40-100
3	P 40	16-40
4	P 16	10-16
5	P 1.6	1.0-1.6