


DURAN® Gas Distribution Tube




Typical application: reactions between gases and liquids.

Cat. No.	d (OD) (mm)	d ₁ (OD) (mm)	h (mm)	Porosity	Pack Unit
side mounted filter cup					
258560002	6	22	250	0	5
258560105	6	22	250	1	5
258560208	6	22	250	2	5
centrally mounted filter cup					
258561007	11	34	250	1	5
258561101	9	25	250	1	5
258561204	9	25	250	2	5
258562106	11	34	250	1	5



DURAN® Gas Washing Bottle

Drechsel type head



With screw-connection system. The insertion height of the head is adjustable. Individual parts can also be ordered separately.

Typical applications: cleaning ("washing") of gases with solvents.


Cat. No.	Capacity (mL)	DIN Thread (GL)	Porosity	Hose connection d (OD) (mm)	Bowl (mm)	Pack Unit
without filter disk						
247130008	500	45		9		1
with filter disk						
257040101	500	45	1	9	25	1



Cat. No.	Description	Pack Unit
Individual parts		
247130205	Drechsel-head, without filter disc	1
257540109	Drechsel-head, with filter disc (Por. 1)	1
218014401	DURAN® Laboratory bottle 500 mL, with DIN thread GL 45	10
292550603	Plastic hose connection, straight	10
292270508	Screw-caps with aperture, PBT, GL 14	10
292271007	Screw-caps with aperture, PBT, GL 45	10
292282501	VQM rubber ring (26 x 42 x 5 mm)	10

DURAN® Gas Washing Bottle

head with filter disk, with standard ground joint



Typical applications: cleaning ("washing") of gases with solvents.

Cat. No.	Description	Capacity (mL)	h (mm)	h ₁ (mm)	Neck	Porosity	Hose connection d (OD) (mm)	Bowl (mm)	Pack Unit
257010107		100	250	180	34/35	1	10	25	1
257020108		250	250	160	45/40	1	10	34	1
Individual parts									
257520107	Gas washing bottle, head with fritted disc (100 mL)					1			1
257520108	Gas washing bottle, head with fritted disc (250 mL)					1			1

