

With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. With proven DURAN® properties. Complete with blue screw cap (PP, integral lip seal) and pouring ring (PP) for drip-free pouring and clean, safe working. Service temperature level of screw cap and pouring ring: +140 °C.

Typical applications: storage, sample preparation, transport. Autoclaving media.

Cat. No.	Capacity (mL)	DIN Thread (GL)	d (OD) (mm)	h (mm)	Remark	Pack Unit
with screw cap and pouring ring from PP (blue)						
218010851	10	25	36	54	Acceptance within ISO 4796 standard has been requested. With specially shaped glass edge for improved pouring out (so that an additional plastic pouring ring is not required).	10
218011453	25	25	36	74	With specially shaped glass edge for improved pouring out (so that an additional plastic pouring ring is not required).	10
218011753	50	32	46	91		10
218012458	100	45	56	105		10
218012955	150	45	62	115		10
218013651	250	45	70	143		10
218014459	500	45	86	181		10
218015155	750	45	95	208		10
218015455	1 000	45	101	230		10
218016357	2 000	45	136	265		10
218016957	3 500	45	160	300		1
218017353	5 000	45	182	335		1
218018658	10 000	45	227	415		1
218018855	15 000	45	268	450		1
218019157	20 000	45	288	510		1
218019251	25 000	45	316	545		1

DURAN® Original Laboratory Bottle

with DIN thread

