


DURAN® Protect Laboratory Bottle

with DIN thread, plastic coated



With easy-to-read scale. In fired-on, highly durable white ceramic. Service temperature limit of the PU plastic coating: -30°C to $+135^{\circ}\text{C}$. The coating provides scratch, leak* and splinter* protection and is ideally suited to both the transport and storage of hazardous media or valuable samples. UV protection up to approx. 380 nm wavelength. High transparency. Suitable for microwaving. (* only applies to bottles 5 000 mL and less)

Typical applications: storage, transport and safe handling of hazardous or valuable substances.

Cat. No.	Capacity (mL)	DIN Thread (GL)	d (OD) (mm)	h (mm)	Remark	Pack Unit
with screw cap and pouring ring						
218052453	100	45	56	100		10
218052959	150	45	62	110		10
218053655	250	45	70	138		10
218054454	500	45	86	176		10
218055159	750	45	95	203		10
218055459	1 000	45	101	225		10
218056352	2 000	45	136	260		10
218056952	3 500	45	160	295		1
218057357	5 000	45	182	330		1
without screw cap and pouring ring						
218050806	10	25	36	50	Acceptance within ISO 4796 standard has been requested. With specially shaped glass edge for improved pouring out (so that an additional plastic pouring ring is not required).	10
1092676	25	25	36	70	With specially shaped glass edge for improved pouring out (so that an additional plastic pouring ring is not required).	10
1092677	50	32	46	87		10
218052404	100	45	56	100		10
218052901	150	45	62	110		10
218053606	250	45	70	138		10
218054405	500	45	86	176		10
218055101	750	45	95	203		10
218055401	1 000	45	101	225		10
218056303	2 000	45	136	260		10
218056903	3 500	45	160	295		1
218057308	5 000	45	182	330		1
218058604	10 000	45	228	410		1
218058801	15 000	45	268	445		1
218059103	20 000	45	289	505		1