



**Reactor Holder Type Schmizo, Polyester
(Glass-Fiber Reinforced)**

SKU: 12.344.37

Product Description

Max. Temperature Resistance 160°C (temporary up to 200°C)

Reactor Holders DN60, DN100 and DN120:

- incl. thread rod M8 13 x 200 mm (see picture 1)
- can be modified with a twin rod (see picture 3)

Reactor Holder DN 150:

- with two lateral support segments (see picture 2)
- for rods up to 13 mm diameter (not included)

Variations

SKU	DN	Flanschttyp
12.344.50	120	SCHMIZO
12.344.80	150	SCHMIZO
12.344.98	60	SCHMIZO
12.346.37	100	SCHOTT
12.346.50	120	SCHOTT
12.346.80	150	SCHOTT
12.346.98	60	SCHOTT
12.344.37	100	SCHMIZO

Customization

We would be happy to discuss with you the development of customer-specific functions, to meet your specific process requirements. Thanks to our technical knowledge and skills in manufacturing scientific glass devices, we can implement almost any requirement and ensure that the customized reactor features are both reliable and user-friendly.