

Accessories for flat flange vessels.

| Cat. No.  | DN  | Pack Unit |
|-----------|-----|-----------|
| 290713407 | 60  | 1         |
| 290714609 | 100 | 1         |
| 290715108 | 120 | 1         |
| 290715708 | 150 | 1         |
| 290716104 | 200 | 1         |

## Quick Release Clamp

from stainless steel, with retaining clip




For secure fitting of the lid or the reaction vessel using two clamping rods.

| Cat. No.  | DN  | Pack Unit |
|-----------|-----|-----------|
| 290734602 | 100 | 1         |
| 290735701 | 150 | 1         |

## Holding Device for Reaction Vessels

from chrome-nickel steel



| Cat. No.  | d (OD) (mm) | h (mm) | Cone (NS) | Pack Unit |
|-----------|-------------|--------|-----------|-----------|
| 215704207 | 40          | 240    | 29/32     | 1         |
| 215704807 | 54          | 315    | 34/35     | 1         |

## DURAN® Calcium Chloride Cylinder

