

With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable ceramic. Complete system comprising: DURAN® pressure plus+ laboratory bottle with 4-port screw cap (PP); four screw connections (black, M8 thread) and silicone seals. Connection of different hose diameters (1.6 mm and 3.2 mm) as well as sterile pressure equalisation sets (syringe filter 0.2 µm) is possible. Unused ports can be sealed with silicone blanking seals.

Typical applications: safe transfer of liquid media within a closed and sterile system (evaporation is reduced).

Cat. No.	Description	Capacity (mL)	DIN Thread (GL)	d (OD) (mm)	h (mm)	Pack Unit
1129821	HPLC Bottle complete	500	45	86	176	2
1129820	HPLC Bottle complete	1 000	45	101	225	2
Accessories						
1129812	Screw cap HPLC, GL 45, 4 ports, complete (GL 45 screw cap, 4 x M8 black screw caps, 12 x silicone seals for 1.6, 3.2 mm tubing or blanks)					2
1129813	Replacement parts for HPLC screw cap includes M8 caps and 1.6 , 3.2 mm and blind silicone seals					1
1137801	Pressure compensation set 4-port cap (incl. 0.2 µm syringe filter)					1
1129819	Spare syringe filter for pressure compensation, 0.2 µm					2

## DURAN® HPLC Bottle

with DIN thread, GL 45

