


DURAN® Original GL Screw Cap from PP

with lip seal



A
121 °C

Tmax.
140 °C

Available in the colours blue, green, yellow and grey with matching colour pouring rings. Distinguishing media types is simplified and interchanging of screw caps and carry over of substances is practically eliminated.

Typical applications: visual aid to the safe identification of different media.

Cat. No.	DIN Thread (GL)	d (OD) (mm)	h (mm)	Colour	Pack Unit
screw cap					
292391307	25	33	19	blue	10
292391907	32	40	24	blue	10
292392809	45	54	25	blue	10
293382802	45	54	25	yellow	10
293382868	45	54	25	green	10
293382884	45	54	25	grey	10
Pouring ring					
292421907	32		4	blue	10
292422809	45		4	blue	10
1089911	45		4	green	10
1089914	45		4	grey	10
1089917	45		4	yellow	10

ס.י.מדה בע"מ
יבוא ושיווק ציוד ומכשור למעבדות

