

High leak tightness through use of PTFE coated silicone seal (peroxide-cured silicone). More chemically and heat resistant than PP screw cap. A matching ETFE pouring ring is also available, permitting clean, drip-free use.


Cat. No.	DIN Thread (GL)	d (OD) (mm)	h (mm)	Pack Unit
screw cap				
292400806	14	20	17	10
292401108	18	23	20	10
292401305	25	33	23	10
292401905	32	41	26	10
292402807	45	54	28	10
Pouring ring				
292441909	32		4	10
292442802	45		4	10

ס.י.מדה בע"מ
יבוא ושיווק ציוד ומכשור למעבדות



DURAN® Red GL High Temperature Screw Cap from PBT

with PTFE coated silicone seal



A
121 °C

Tmax.
180 °C