


For 10 microscope slides 76 x 26 mm. Note: Do not clean staining dishes and staining jars at temperatures above 60 °C (glass corrosion is possible).

Cat. No.	h (mm)	l (mm)	b (mm)	Pack Unit
233140001	90	90	40	10

Staining Dish Hellendahl Type

from Soda-lime Glass, straight sided




For 16 microscope slides 76 x 26 mm, with widening towards the top. Note: Do not clean staining dishes and staining jars at temperatures above 60 °C (glass corrosion is possible).

Cat. No.	h (mm)	l (mm)	b (mm)	Pack Unit
233150002	100	60	60	10

Staining Dish Hellendahl Type

from Soda-lime Glass




For 10 microscope slides 76 x 26 mm. Note: Do not clean staining dishes and staining jars at temperatures above 60 °C (glass corrosion is possible).

Cat. No.	h (mm)	l (mm)	b (mm)	Pack Unit
233160003	40	90	70	10

Staining Dish Schiefferdecker Type

from Soda-lime Glass



Glass Box from Soda-lime Glass

for staining tray 213170003



Note: Do not clean staining dishes and staining jars at temperatures above 60 °C (glass corrosion is possible).

Cat. No.	h (mm)	l (mm)	b (mm)	Pack Unit
233180005	70	108	90	10

DURAN® Staining Tray

For 10 microscope slides 76 x 26 mm or each width up to 52 mm.



Cat. No.	h (mm)	l (mm)	b (mm)	Pack Unit
213170003	70	88	40	10



Stainless Steel Handle

for staining tray

Cat. No.	Pack Unit
290750002	10

