


Cat. No.	d (OD) (mm)	h (mm)	Volume approx. (mL)	Wall thickness (mm)	Pack Unit
231720184	9.75	75	4	0.8	814
231720365	10.00	75	4	0.6	766
231720587	11.75	75	5	0.8	550
231720862	11.75	75	6	0.55	550
231720965	12.25	75	7	0.55	500
231720998	12.25	75	6	0.8	500
231720784	11.75	100	8	0.8	550
231721197	12.25	100	9	0.8	500
231721489	15.50	100	14	0.8	310
231721094	15.75	100	15	0.9	310
231721283	12.25	120	18	0.8	500
231721883	10.00	150	8	0.8	766
231721986	15.50	150	19	0.8	310
231722193	15.50	160	22	0.8	310

Disposable Culture Tube

from Soda-lime Glass, straight rim



ס.י.מדה בע"מ

יבוא ושיווק ציוד ומכשור למעבדות

WWW.SIMADA.CO.IL 03-5591008



DWK
LIFE SCIENCES