



## Vials

Clear and amber dram sample vials are made of Type I borosilicate glass, providing exceptional resistance to heat shock and chemical leaching.

Finneran dram vials are available

precleaned and precleaned/ certified. Please contact us for ordering information.

### Case Pack - 100 pieces

Cat. No.	Description
34013-1545	1 Dram, 15x45mm Clear Vial, 13-425mm Thread (1,000 per case)
34013-1545A	1 Dram, 15x45mm Amber Vial, 13-425mm Thread (1,000 per case)
38015-1760	2 Dram, 17x60mm Clear Vial, 15-425mm Thread
38015-1760A	2 Dram, 17x60mm Amber Vial, 15-425mm Thread
312015-1965	3 Dram, 19x65mm Clear Vial, 15-425mm Thread
312015-1965A	3 Dram, 19x65mm Amber Vial, 15-425mm Thread
316018-2170	4 Dram, 21x70mm Clear Vial, 18-400mm Thread
324020-2385	6 Dram, 23x85mm Clear Vial, 20-400mm Thread

## Screw Thread Closures

### Shelf Pack - 100 pieces; Case Pack - 1000 pieces

Cat. No.	Description
5310-13	13-425mm Black, Polypropylene Open Hole Cap
5310-13W	13-425mm White, Polypropylene Open Hole Cap
5310-15	15-425mm Black, Polypropylene Open Hole Cap
5310-15W	15-425mm White, Polypropylene Open Hole Cap
5310-18	18-400mm Black, Polypropylene Open Hole Cap
5310-20	20-400mm Black, Polypropylene Open Hole Cap
5320-13	13-425mm Solid Top, Black Polypropylene Unlined
5320-13W	13-425mm Solid Top, White Polypropylene Unlined
5320-15	15-425mm Solid Top, Black Polypropylene Unlined
5320-15W	15-425mm Solid Top, White Polypropylene Unlined
5320-18	18-400mm Solid Top, Black Polypropylene Unlined
5320-20	20-400mm Solid Top, Black Polypropylene Unlined
5360-13	13-425mm Solid Top, Black Polypropylene PTFE/F217 Lined
5360-13W	13-425mm Solid Top, White Polypropylene PTFE/F217 Lined
5360-15	15-425mm Solid Top, Black Polypropylene PTFE/F217 Lined
5360-15W	15-425mm Solid Top, White Polypropylene PTFE/F217 Lined

## Closures

### Septa for Screw Thread Closures

**Shelf Pack - 100 pieces; Case Pack - 1000 pieces**

<b>Cat. No.</b>	<b>Description</b>
606050G-13	13mm x 0.060" Gray PTFE/SIL Septa for Dram Vials
606050G-15	15mm x 0.060" Gray PTFE/SIL Septa for Dram Vials
606050G-18	18mm x 0.060" Gray PTFE/SIL Septa for Dram Vials
606050G-20	20mm x 0.060" Gray PTFE/SIL Septa for Dram Vials

### Preassembled Convenience Packs

Dram vials assembled with black polypropylene closures and septa.

**Case Pack - 100 pieces**

<b>Cat. No.</b>	<b>Description</b>
82020-1235	1/2 Dram, 12x35mm Clear Vial, 8-425mm Solid Top Cap, PTFE/F217 Lined
84020-1545	1 Dram, 15x45mm Clear Vial, 13-425mm Solid Top Cap, PTFE/F217 Lined
88020-1760	2 Dram, 17x60mm Clear Vial, 15-425mm Solid Top Cap, PTFE/F217 Lined
812020-1965	3 Dram, 19x65mm Clear Vial, 15-425mm Solid Top Cap, PTFE/F217 Lined
812020-1965A	3 Dram, 19x65mm Amber Vial, 15-425mm Solid Top Cap, PTFE/F217 Lined
812030-1965	3 Dram, 19x65mm Clear Vial, 15-425mm Open Hole Cap, PTFE/SIL Lined
812030-1965A	3 Dram, 19x65mm Amber Vial, 15-425mm Open Hole Cap, PTFE/SIL Lined
816020-2170	4 Dram, 21x70mm Clear Vial, 18-400mm Solid Top Cap, PTFE/F217 Lined
816030-2170	4 Dram, 21x70mm Clear Vial, 18-400mm Open Hole Cap, PTFE/SIL Lined
824020-2385	6 Dram, 23x85mm Clear Vial, 20-400mm Solid Top Cap, PTFEF217 Lined