

TRIDENT VIALS



Screw Neck Vials

Neutral glass, tall form, without closure. Polypropylene caps must be ordered separately

Code	Nominal capacity dram	ml	Height mm	Diameter mm	Neck bore mm	Quantity per pack	Trident code	Cap reference
T101/V1	0.5	1.75	36	11.0	5.5	882	09/001	T001/C1
T101/V2	1.0	3.50	46	12.5	5.5	666	09/002	T001/C1
T101/V3	2.0	7.0	58	17.0	9.0	399	09/009	T001/C3
T101/V4	3.0	10.5	67	17.0	9.0	399	09/008	T001/C3
T101/V5	4.0	14.0	72	20.0	13.0	264	09/007	T001/C5
T101/V6	6.0	21.25	86	21.0	13.0	231	09/006	T001/C5
T101/V7	8.0	28.25	96	23.0	15.5	190	09/005	T001/C7

Screw Neck Vials

Neutral glass, squat form, without closure. Polypropylene caps must be ordered separately

Code	Nominal capacity dram	ml	Height mm	Diameter mm	Neck bore mm	Quantity per pack	Trident code	Cap reference
T102/V1	2.0	7.0	42	20.0	13.0	264	09/003	T001/C5
T102/V2	4.0	14.0	58	23.0	15.5	190	09/004	T001/C7
T102/V3	8.0	28.25	72	27.5	21.0	140	09/010	T002/C3

Screw Caps

Black polypropylene with foil-faced wads

Code	Type
T001/C1	To fit T101/V1 & T101/V2
T001/C3	To fit T101/V3 & T101/V4
T001/C5	To fit T101/V5, T101/V6 & T102/V1
T001/C7	To fit T101/V7, & T102/V2
T002/C3	To fit T102/V3



Rolled Rim Vials

Neutral glass, without closure

These vials will accept either a polythene push-in plug or a clip-on closure which must be ordered separately

Code	Nominal capacity dram	ml	Height mm	Diameter mm	Neck bore mm	Quantity per pack	Trident code
T103/V1	2.0	7.0	34	23.0	17.15	190	07/014
T103/V2	3.0	10.5	46	23.0	17.15	190	07/013
T103/V3	4.0	14.0	48	25.5	17.15	162	07/010
T103/V4	6.0	21.25	66	25.5	17.15	162	07/011
T103/V5	8.0	28.25	72	27.5	17.15	140	07/012

Caps for Rolled Rim Vials

Code	Type
T003/C1	Clip-on closure, polythene
T004/C1	Push-in plug, polythene



Wide Mouth Vials

Neutral glass, tall form, without closure

Polythene plugs must be ordered separately

Caps for Wide Mouth Vials

Code	Type
T005/C1	Push-in plug, polythene to fit T105/V1 & V2
T005/C3	Push-in plug, polythene to fit T105/V3

Code	Nominal capacity dram	ml	Height mm	Diameter mm	Neck bore mm	Quantity per pack	Trident code	Cap reference
T105/V1	0.5	1.75	36	11.0	8.25	882	08/016	T005/C1
T105/V2	1.0	3.5	46	12.5	8.25	666	08/017	T005/C1
T105/V3	2.0	7.0	58	17.0	10.25	399	08/018	T005/C3