



Magnetic Stirrer - MICROSTIRRER

GENERAL FEATURES	
Construction material:	epoxy painted metal structure
Power:	0.6 W
Weight:	0.55 Kg (1.1 lb)
Dimensions (WxHxD):	120x48x128 mm (4.7x1.9x5 in)
PERFORMANCE	
Electronic speed regulation:	up to 1100 rpm
Stirring volume (H ₂ O):	up to 5 liters
Speed control:	excellent speed control even at low revolutions

High Reliability

- Control panel is shielded from liquid spills
- Epoxy-coated tough steel casing
- Long life operation
- Accepts up to **5 liter** flasks

Applications in Mind

- Pale surface for visualization of particles and color changes
- **Remains cold**, even after days of continuous use
- Specifically designed for microtitration
- Excellent speed control, even at low speed

Simple & Practical

- Speed regulation up to **1100 rpm**
- Universal mains voltage / frequency range
- Light weight