



## Hot Plate - RC Series

**Single-position and two-position heating plates with aluminum top and temperature regulation.**

### Excellent Performance

- **Uniform** heat distribution with coated aluminum alloy heating plate
- Excellent resistance to chemicals
- Maximum temperature **370 °C**
- Accepts up to **15 liter** flasks (Ø 155 mm hotplate)

### Simple & Practical

- Control panel is IP 42 protected from liquid spills with run-off groove (RC only)
- Choice of mains voltage / frequency model





### Ergonomic Design

- Inclined control panel for easy access and visibility (RC)
- Control panel shielded against spillage (RC 2)
- **AluBlocks™** maximize temperature control & save handling time (RC only)
- Run different tubes by combining **AluBlocks™** (RC only)

GENERAL FEATURES	
Construction material:	epoxy painted aluminum structure
Heating plate:	aluminum alloy coated with a special protection
Heating plate diameter:	155 mm
Models:	<b>RC</b> single-position <b>RC2</b> two-position

Power:	<b>RC</b> 600 W <b>RC2</b> 1200 W
Weight:	<b>RC</b> 1.4 Kg (3.1 lb) <b>RC2</b> 3.3 Kg (7.3 lb)
Dimensions (WxHxD):	<b>RC</b> 165x115x280 mm (6.5x4.5x11 in) <b>RC2</b> 340x90x246 mm (13.4x3.5x9.7 in)
<b>PERFORMANCE</b>	
Temperature regulation:	from room temp. to 370 °C

## Product Code

Description	Power Supply	Code
 <b>RC</b>	230 V / 50-60 Hz	F20700431
 <b>RC</b>	115 V / 50-60 Hz	F20710431
 <b>RC2</b>	230 V / 50-60 Hz	F20700430
 <b>RC2</b>	115 V / 50-60 Hz	F20710430

## Accessories



**A00001069**Support rod  
(RC only)



**A00000273**Hemispheric  
bowl for 50 ml flasks



**A00000258**Hemispheric  
bowl for 100 ml flasks



**A00001071**Hemispheric  
bowl for 250 ml flasks



**A00001072**Hemispheric  
bowl for 500 ml flasks



**A00001073**Hemispheric  
bowl for 1000 ml flasks



**A00000228**AluBlock™  
Base (RC only)



**A00000229**Orange  
AluBlock™, 4 pos. Ø28  
x h 43mm (RC only)



**A00000230**Green  
AluBlock™, 4 pos. Ø28  
x h 30mm (RC only)

**A00000231**Black  
AluBlock™, 4 pos. Ø28  
x h 24mm (RC only)

**A00000232**Red  
AluBlock™, 4 pos.  
Ø21,6 x h 31,7 mm (RC  
only)

**A00000233**Blue  
AluBlock™, 8 pos.  
Ø17,8 x h 26mm (RC  
only)

**A00000234**Gold  
AluBlock™, 11 pos.  
Ø15,2 x h 20mm (RC  
only)