

MH-4

Analog Aluminium Top Hotplate Magnetic Stirrer



Easy to clean chemical resistant top plate.

Temperature selection is provided by a Hydraulic Thermostat and the stirring by quiet direct drive motor / magnet system.

Temperature ranges between 50-300 C (122 F- 572 F) and the speed can be set between 100-1000 rpm. Top plate: 17.5 x 17.5 cm. Overall Dim. (W x H x D) 20 x 14 x 26 cm (D includes rear bracket for support rod). Nett weight : 3.2 Kgs. 230 Volt (or 110V), 50/60 Hz , 750 Watt