



Beakers stainless steel

Type 1: with rim

Type 2: with rim and spout

Type 3: with rim and handle

Type 4: with rim, spout and handle

Material	C (Height)	D/Ø	ml	Type	Article no.
18/10 Stainless	60	50	100	1	8500
18/10 Stainless	85	65	250	1	8501
18/10 Stainless	110	80	500	1	8502
18/10 Stainless	130	100	1000	1	8503
18/10 Stainless	180	120	2000	1	8504
18/10 Stainless	180	150	3000	1	8505
18/10 Stainless	210	180	5000	1	8506
18/10 Stainless	60	50	100	2	8510
18/10 Stainless	85	65	250	2	8511
18/10 Stainless	110	80	500	2	8512
18/10 Stainless	130	100	1000	2	8513
18/10 Stainless	180	120	2000	2	8514
18/10 Stainless	180	150	3000	2	8515
18/10 Stainless	210	180	5000	2	8516
18/10 Stainless	60	50	100	3	8520
18/10 Stainless	85	65	250	3	8521
18/10 Stainless	110	80	500	3	8522
18/10 Stainless	130	100	1000	3	8523
18/10 Stainless	180	120	2000	3	8524
18/10 Stainless	180	150	3000	3	8525
18/10 Stainless	210	180	5000	3	8526
18/10 Stainless	60	50	100	4	8530
18/10 Stainless	85	65	250	4	8531
18/10 Stainless	110	80	500	4	8532
18/10 Stainless	130	100	1000	4	8533
18/10 Stainless	180	120	2000	4	8534
18/10 Stainless	180	150	3000	4	8535
18/10 Stainless	210	180	5000	4	8536