



- Propellerrührer, 2 Flügel flexibel
- Propeller stirrers, 2 flexible blades
- Tiges à hélice, 2 pales basculantes
- Hélice agitadora, 2 palas móviles




Material	A	D	D1	Artikel Nr.
18/10-STAHL	300	8	60	9603
18/10-STAHL	400	8	60	9604
18/10-STAHL	500	8	60	9605

- Propellerrührer, 2 Flügel fix
- Propeller stirrers, 2 fix blades
- Tiges à hélice, 2 pales fixes
- Hélice agitadora, 2 palas fijas




Material	A	D	D1	Artikel Nr.
18/10-STAHL	300	8	50	9703
18/10-STAHL	400	8	50	9704
18/10-STAHL	500	8	50	9705

- Ankerrührer
- Anchor stirrers
- Tiges à ancre
- Hélice forma ancla





Material	A	D	D1	Artikel Nr.
18/10-STAHL	500	8	70	9610
18/10-STAHL	650	10	90	9611
18/10-STAHL	800	10	140	9612

- Propellerrührer, 3 Flügel
- Propeller stirrers, 3 blades
- Tiges à hélice, 3 pales
- Hélice agitadora, 3 palas




Material	A	D	D1	Artikel Nr.
18/10-STAHL	300	8	50	9103
18/10-STAHL	400	8	50	9104
18/10-STAHL	500	8	50	9105
18/10-STAHL	300	8	70	9109
18/10-STAHL	400	8	70	9110
18/10-STAHL	500	8	70	9111
18/10-STAHL	300	10	100	9112
18/10-STAHL	400	10	100	9113
18/10-STAHL	500	10	100	9114
18/10-STAHL	650	10	70	9115
18/10-STAHL	800	10	100	9116

- Propellerrührer, 4 Flügel vertikal
- Propeller stirrers, 4 vertical blades
- Tiges à hélice, 4 pales verticales
- Hélice agitadora, 4 palas verticales



Material	A	D	D1	Artikel Nr.
18/10-STAHL	300	8	50	9053
18/10-STAHL	400	8	50	9054
18/10-STAHL	500	8	50	9055
18/10-STAHL	300	10	100	9056
18/10-STAHL	400	10	100	9057
18/10-STAHL	500	10	100	9058

- Propellerrührer, 4 Flügel
- Propeller stirrers, 4 blades
- Tiges à hélice, 4 pales
- Hélice agitadora, 4 palas



Material	A	D	D1	Artikel Nr.
18/10-STAHL	300	8	50	9003
18/10-STAHL	400	8	50	9004
18/10-STAHL	500	8	50	9005
18/10-STAHL	300	8	100	9009
18/10-STAHL	400	8	100	9010
18/10-STAHL	500	8	100	9011
18/10-STAHL	500	8	70	9012
18/10-STAHL	650	10	100	9013
18/10-STAHL	800	10	100	9014

### ■ Zentrifugalrührer, 2 flexible Flügel

für Gefässe mit engem Hals

### ■ Centrifugal stirrers, 2 flexible blades


for jars with tight neck

### ■ Tiges à hélice centrifuge, 2 pales flexibles

pour récipients avec col étroit



### ■ Hélice centrifuga 2 palas flexibles

para recipientes con cuello estrecho



Material	A	D	D1	Artikel Nr.
18/10-STAHL	300	8	90/15	9209
18/10-STAHL	400	8	90/15	9210
18/10-STAHL	500	8	90/15	9211
18/10-STAHL	300	10	90/15	9212
18/10-STAHL	400	10	90/15	9213
18/10-STAHL	500	10	90/15	9214
18/10-STAHL	650	10	90/15	9215
TEFLON COATING	500	8	90/15	9216

- Flächenrührer, 3 Loch
- Paddle stirrers, 3 holes
- Tiges à palette, 3 trous
- Hélice de pala, 3 agujeros




Material	A	D	D1	Artikel Nr.
18/10-STAHL	300	8	70	9403
18/10-STAHL	400	8	70	9404
18/10-STAHL	500	8	70	9405
18/10-STAHL	300	10	70	9406
18/10-STAHL	400	10	70	9407
18/10-STAHL	500	10	70	9408

### ■ Flächenrührer, 6 Loch

### ■ Paddle stirrers, 6 holes

### ■ Tiges à palette, 6 trous

### ■ Hélice de pala, 6 agujeros



Material	A	D	D1	Artikel Nr.
18/10-STAHL	300	8	70	9503
18/10-STAHL	400	8	70	9504
18/10-STAHL	500	8	70	9505
18/10-STAHL	300	10	70	9506
18/10-STAHL	400	10	70	9507
18/10-STAHL	500	10	70	9508