

## Tubes with Screw Cap



### Polystyrene Tubes


#### Screw Cap, Conical / Round Bottom

- ↳ Further Closures p. 5 | 13
- ↳ Overview: See Technical Appendix for max. relative centrifugal force (RCF)
- ↳ Cell Culture Tubes p. 1 | 21
- ↳ Tube No. 163 177 also available as polypropylene version with white screw cap (Cat.-No. 163 270)



Cat.-No.	164 180	164 161	163 177
Bottom design	conical	conical	round
ø [mm] x height [mm]	16.8 x 100	16.8 x 100	17 x 100
Volume [ml]	12	12	12
Sterile	-	+	-
Screw cap	-	-	black
Bayonet cap <sup>*)</sup>	white	blue	-
Graduation	-	-	-
Writing area	-	-	-
Support skirt	+	+	-
Quantity per bag/case	1300	25/1000	900

<sup>\*)</sup> open by a 1/3 turn



Cat.-No.	188 161	188 171	186 161	186 171
Bottom design	conical	conical	round	round
ø [mm] x height [mm]	17 x 120	17 x 120	17 x 120	17 x 120
Volume [ml]	15	15	15	15
Sterile	+	+	+	+
Screw cap	blue	blue	blue	blue
Graduation	+	+	-	-
Writing area	+	+	-	-
Support skirt	-	-	-	-
Quantity per bag/case	50 <sup>*)</sup> /500	100/1000	50 <sup>*)</sup> /500	100/1000

<sup>\*)</sup> box

<sup>\*)</sup> box

1 Cell/ Tissue Culture  
2 HTS- Microplates  
3 Immunology/ HLA  
4 Microbiology/ Bacteriology  
5 Tubes/Multi- Purpose Beakers  
6 Liquid Handling  
7 Molecular Biology  
8 Protein Crystallisation  
9 Separation  
10 Biochips/ Microfluidics  
11 Cryo-Technics  
12 Lids/Sealers/ CapMats  
13 Reaction Tubes/ Analyser Cups  
14 Accessories