



Specimen Containers, 30ml , Universal, Screw Cap

1. Manufactured from virgin polystyrene and polypropylene. All materials are non-cytotoxic
2. CE Marked in accordance with the European Directive 98/79/EC
3. For In Vitro Use only
4. Leakproof-tested in accordance with BS5213:1975(1989) Appendix A
5. Unique leak tight cap ensures exceptional leak proof seal
6. Available pre-filled with boric acid (0.4g) for the preservation of urine samples
7. Containers are available in a full range of label options, no label, plain label and printed label (for the clinical market)

C/N	Desc. of Item	Sterility	Inner packing	Case Qty.
22100	PS, No label	EO	50/zip bag	400
22101	PS, Printed label	EO	50/zip bag	400
22102	PS, Plain label	EO	50/zip bag	400
22103	PS, with Boric Acid, Printed label	EO	50/zip bag	400
22110	PP, No label	EO	50/zip bag	400
22111	PP, Printed label	EO	50/zip bag	400
22112	PP, Plain label	EO	50/zip bag	400
22113	PP, with Boric Acid, Printed label	EO	50/zip bag	400